BRAUN

Pulse oximeter 1


YK-81CEU

GB	English.....	2
DE	Deutsch	12
DK	Dansk.....	23
ES	Español	33
FI	Suomi	43
FR	Français	53
GR	Ελληνικά	63
IT	Italiano.....	73
NL	Nederlands.....	83
NO	Norsk.....	93
PL	Polski.....	103
PT	Portugues	113
SA	اللغة العربية.....	131
SE	Svenska.....	133
TR	Türkçe.....	143

Product description

OLED display modes

Not made with natural rubber latex

Naturkautschuk

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naturgummilætex

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Thank you for purchasing the Braun Pulse oximeter, YK-81CEU. This pulse oximeter is a high-quality product tested in accordance with international ISO safety and performance standards. This device uses red and infrared light at set wavelengths along with a photodetector to measure oxygen saturation in blood and the associated pulse rate.

The Braun Pulse oximeter is intended to be used by consumers in a home-use setting (such as home setting or places of business).

Please read these instructions carefully before using this product and keep the instructions and the pulse oximeter in a safe place.

Indications for use

The Braun Pulse oximeter (YK-81CEU) is a non-sterile, reusable device intended for spot checking or intermittent determination of functional oxygen saturation in arterial blood, as well as the associated pulse rate by using the tip of the finger as the measurement site. The device is only intended for people 18 years of age or older, in a home-use environment (such as a home setting or places of business).

The Pulse oximeter displays the %SpO₂, pulse rate, as well as the pulse amplitude using a vertical bar graph.

The Pulse oximeter should only be used by people 18 years of age or older.

Warnings and precautions

Please read the manual carefully before use.

Device isn't suitable for continuous monitoring.

The measurement results of Braun Pulse oximeter cannot replace the doctor's diagnosis, please follow the guidance of a professional doctor.

Pulse oximeters are sensitive to motion. Keep hands still while taking a reading.

Pulse oximeters require sufficient blood flow to obtain proper readings. Poor blood circulation can result in inaccurate readings. If your hands are cold or if you have poor circulation, warm your hands by rubbing them together or use another warming method before attempting to obtain a reading. A tourniquet, blood pressure cuff, other devices or clothing that cause blood flow restrictions may also result in inaccurate readings.

Fingernail polish or acrylic nails obstruct the light transmission through the fingertip and may also result in inaccurate readings.

Your finger must be clean for a proper reading.

The Pulse oximeter should be clean in order to get a proper reading. Please ensure the finger clip area is cleaned, following the cleaning instruction in the manual.

If a reading is difficult to obtain, switch to another finger or to the other hand. The index or middle finger is recommended for taking measurements.

There are a number of other conditions which may lead to an inaccurate reading; including but not limited to, recent medical tests that included an injection of dyes, use of arterial catheters, a weak pulse, low levels of hemoglobin in the blood, low perfusion (the quality of your pulse), elevated levels of dysfunctional hemoglobin, the strength and type of light that you are in while using the pulse oximeter, and the existence of cell phones, radios, and fixed transmitters within certain ranges of the pulse oximeter during use.

This device does not provide a low SpO₂ alarm.

DO NOT use on the same finger continuously for over 2 hours.

The device is not appropriate for the constant monitoring of the user. Prolonged and continuous monitoring may increase the potential of unexpected changes of dermal conditions such as putrescence, etc. A feeling of discomfort or pain may occur if using the device continuously, especially for those users who experience poor perfusion or immature dermatographia by light collimation.

In addition to the items described in the Warnings and precautions section, inaccurate measurements may be caused by factors including but not limited to:

- Immersing the device in liquid.
- Significant levels of dysfunctional hemoglobin (such as carboxyhemoglobin or methemoglobin).
- Intravascular dyes such as indocyanine green or methylene blue.
- Measurements may be adversely affected in the presence of high ambient light such as direct sunlight. In bright light conditions, cover the sensor area if necessary.
- Excessive user movement.
- High-frequency electrosurgical interference and defibrillators.
- Venous pulsations.
- Hypotension, severe vasoconstriction, severe anemia, or hypothermia.

DO NOT modify this equipment without the authorization of the manufacturer.

Fingertip pulse oximeter

General description

Oxygen mixes with blood in the lungs and is transported in blood mainly by hemoglobin. Chemically the oxygen molecule combines loosely and reversibly with the heme portion of the hemoglobin. The amount of oxygen that combines with the hemoglobin is determined by the partial pressure of oxygen. The amount of oxygen in blood, expressed as a percentage of hemoglobin saturation is a measure of oxygen (SaO₂) in the artery and is estimated by a pulse oximeter as an SpO₂ reading.

The amount of oxygen in blood (as measured by the pulse oximeter SpO₂ reading) is a very important physiological parameter in respiratory and circulatory systems. Many respiratory diseases reduce the oxygen-hemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation

Malfunction caused by anesthesia, trauma resulting from major surgery or medical examination, or other diseases due to virus or bacteria can also cause a reduction in oxygen hemoglobin saturation in arterial blood which can be detected by the pulse oximeter through a low SpO₂ reading. This can result in symptoms such as depression, vomiting, and asthenia in patients.

The fingertip pulse oximeter is a compact, low-power consumption, portable device. When the index finger is inserted into the fingertip photoelectric sensor, the display screen of the Braun Pulse oximeter will show the measured value of %SpO₂ and pulse rate as well as a vertical bar graph showing the pulse amplitude.

Measurement principle

Pulse oximetry is based on the difference in the selective absorption of red light by oxygenated blood or oxyhemoglobin (HbO₂, which is bright red) as compared to deoxyhemoglobin (Hb, which is dull red) or de-oxygenated blood. This selective absorption is based on the Beer Lambert Law. A red (660nm) and infrared LED (940nm) on the top portion of the prospective clamp finger-type sensor emits light at these two frequencies. The red light transmitted through the finger is selectively absorbed by the oxygenated arterial blood and the remaining light is detected by a photodetector on the lower portion of the clamp finger. The infrared light then passes through the finger without any selective absorption and is detected by the photodetector on the lower side of the clamp. The transmitted light captured by the photodetector is converted into a voltage or current signal and then passed through signal conditioning and filtering circuits using a microprocessor. %SpO₂ is then derived by analyzing the changes in the amplitude of these pulsatile signals derived from the photodetector and calibrating them using pulse scanning and recording technology. The pulse rate is then derived from these associated waveforms.


Product description

1. OLED display
2. Finger clip
3. Power/Menu button
4. Lanyard


How to use your Braun Pulse oximeter

Installing and replacing batteries (initial set up)


Push on the edge of the battery cover behind the arrows to slide it from the device.

Install 2 AAA batteries.


Replace battery cover.


Taking a measurement

Press on top and bottom housing to open rubber finger clip.


Insert finger with nail facing up into rubber finger clip.

It is recommended that you use the oximeter on your index finger.

Make sure your finger is fully inserted into the device and not pressed hard into a surface.

Press Power/Menu button on the top of the device.

Keep finger steady while taking a measurement.


Rotating the display

If you want to change the display direction, press and release (approx. 1 second) the Power/Menu button on the front panel. The display can be viewed in 6 different orientations.


Care and cleaning

Use 70% alcohol to clean the pulse oximeter finger insert area and your finger, before and after each test.

Parameter set up

Press the Power/Menu button to turn on the device. With the device turned on, press and hold the Power/Menu button to change the display mode.

Press and hold the Power/Menu button (approx. > 3seconds) until you enter the parameter setup menu 1 and hold down again (approx. 3 seconds) to enter screen 2.

Settings	
Alm Setup	*
Alm	off
Beep	off
Demo	off
Restore	OK
Brightness	3
Exit	

Settings	
Sounds Setup	*
Spo2 Alm Hi	100
Spo2 Alm Lo	85
PR Alm Hi	130
PR Alm Lo	50
+/-	+
Exit	

1. In menu 1 or 2 press the Power/Menu button to cycle through the various parameters. You will see the * move to each parameter as you cycle through. When you are at the parameter you wish to change, press the Power/Menu button for 3 seconds to change the parameter setting.
2. If you wish to restore to factory settings, move the * to the restore selection on menu 1 and hold the Power/Menu button down for approximately 3 seconds.
3. To exit the menus and return to the measurement screen, move the * to the exit selection on menu 1 or 2 and hold the Power/Menu button for 3 seconds, or simply wait 10-12 seconds without interacting with device and it will return to the measurement screen.

Attention

If using device outdoors or under strong light, please increase the screen brightness to improve legibility.

Choose a lower screen brightness to conserve battery power.

Limited 2 year warranty

A. This 2 year limited warranty applies to repair or replacement of product found to be defective in material or workmanship. This warranty does not apply to damage resulting from commercial, abusive, unreasonable use or supplemental damage. Defects that are the result of normal wear and tear will not be considered manufacturing defects under this warranty. **KAZ IS NOT LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THIS PRODUCT IS LIMITED IN DURATION TO THE DURATION OF THIS WARRANTY.** Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights and you may also have other rights which vary from jurisdiction to jurisdiction. This warranty applies only to the original purchaser of this product from the original date of purchase. **B.** At its option, Kaz will repair or replace this product if it is found to be defective in material or workmanship. **C.** This warranty does not cover damage resulting from any unauthorized attempts to repair or from any use not in accordance with the instruction manual. Refer to page 153 for list of European Contact numbers. Please be sure to specify a model number.

NOTE: IF YOU EXPERIENCE A PROBLEM, PLEASE CONTACT CONSUMER RELATIONS FIRST OR SEE YOUR WARRANTY. DO NOT RETURN THE PRODUCT TO THE ORIGINAL PLACE OF PURCHASE. DO NOT ATTEMPT TO OPEN THE DEVICE HOUSING YOURSELF, DOING SO MAY VOID YOUR WARRANTY AND CAUSE DAMAGE TO THE PRODUCT OR PERSONAL INJURY.

Errors and troubleshooting

Error message	Situation	Solution
SpO2 or PR can not be shown normally	Finger is not fully inserted into device User's oxyhemoglobin value is too low to be measured	Retry by correctly inserting the finger Try a few more times If you can make sure that there is no problem in the product, please go to the hospital in a timely manner for exact diagnosis
SpO2 or PR reading is unsteady	Finger is not fully inserted into device User is moving or finger is shaking	Retry by correctly inserting the finger Remain as still as possible
The oximeter can not be turned on	Inadequate power or power off Batteries are installed incorrectly The oximeter might be damaged	Replace the batteries Please reinstall the batteries Please contact the customer service center
Unit turned off	The product automatically shuts off when no signal can be detected within 8 seconds Inadequate power	Normal Replace the batteries
Low battery displayed	Batteries need to be replaced	Replace the batteries

Product specifications

Type:	Braun Pulse oximeter, YK-81CEU
Display Type:	OLED display
SpO2 Measurement range:	70% - 100%
Accuracy:	80% - 100% ±2% 70% - 79% ±3%
Resolution:	1%
Pulse Rate	
Measurement range:	30BPM~254 BPM
Accuracy:	≤ 100 ±1 BPM > 100 ±2 BPM
Resolution:	1 BPM
Power Supply:	two AAA 1.5V alkaline batteries
Power Consumption:	below 30mAh
Automatically Power-off:	The product automatically shuts off when no signal can be detected within 8 seconds
Dimension:	Approx. 58mm×35mm×30mm
Operating Temperature:	5°C~40°C (41°F ~ 104°F)

Storage Temperature:	-10 °C~40 °C (14 °F ~ 104 °F)
Operating Humidity:	15%~80% RH, Non condensing
Storage Humidity:	10%~95% RH, Non condensing
Atmospheric pressure:	70kPa~106kPa

Maintenance and storage

1. Replace the batteries when low voltage symbol is lit.
2. Clean the surface of the oximeter before it is used to take a measurement.
3. Remove the batteries if you will not use the oximeter for an extended amount of time.
4. For product preservation store the device in -10~40 °C (14-104 °F) and 10%~95% humidity.
5. It is recommended that the product should be kept dry at all times. Excessive humidity might affect its lifetime and even damage the product.
6. This device is factory calibrated and no further calibration is required by the user.

 To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.

 Please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at your local retailer or at appropriate collection points provided in your country

Inquiry

In Europe, any serious incident (e.g., death, life-threatening injury, surgical intervention, etc.) that has occurred in relation to this product should be reported to the manufacturer and/or authorised representative in the European Community and to the competent authority of your member state.

The first 3 numeric digits in the LOT number represent the day of the year of manufacture. The next 2 numeric digits represent the last two numbers of the calendar year of manufacture and the letter(s) at the end designate the manufacturer of the product. (e.g.: LOT No.: 07322YNG this product was made on the day 73, year 2022 at the manufacturer code YNG.)

Box content

- One Braun Pulse oximeter
- One lanyard
- One user manual
- 2 AAA 1.5V alkaline batteries

Symbols and definitions

Type BF Applied Part

Refer to instruction manual/booklet

Manufacturer

3V

Voltage – direct current

Operating Temperature

Storage Temperature

Operating Humidity

Storage Humidity

Importer

Distributor

Batch Code

Serial Number

Medical Device


Catalogue Number

Authorized representative in the European Community

Unique Device Identifier

Atmospheric Pressure

Symbol for compliance with MDD 93/42/EEC requirements

IP22

IP22: Protected against solid foreign objects greater than 12.5 mm in diameter and dripping water when tilted up to 15°

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Xuzhou Yongkang Electronic Science
Technology Co., Ltd
1st & 2nd Floor, 6#01, 6#02, No. 6 Building
1st Phase Economic Development
Manufacturing Zone
LIANDO U Valley, No.6 Leye Road
Xuzhou ETDZ
221000 Xuzhou
PEOPLE'S REPUBLIC OF CHINA

0123

Prolinx GmbH
Brehmstr. 56, 40239
Duesseldorf – Germany

Certain trademarks used under license
from The Procter & Gamble Company or
its affiliates.

© 2022, All rights reserved.

Imported and distributed by:
Kaz Europe Sàrl
A Helen of Troy Company
Q-Center, Chaux 4
CH-1030 Bussigny - Switzerland

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MD **REF** YK-81CEU
P/N: A007999R0
23AUG22

**Helen
of Troy**

**Creative Department
Artwork Specifications**

Brand: Braun
Category: Oximeter
Model: **YK-81CEU**
Artwork Part #: A007999R0
Die Line Part #: NA
Subject: Owners Manual
Region: EMEA
Flat Size: **MM: W160 x H154**

Folded Size: **MM: W80 x H154**

Scale: 1/1
Material: 50# matte text

Page count: 156
Revision: 6
Date: **23AUG22**
Release Date: **24AUG22**
Rerelease Date:

Colors:

Dielines (Do not print)

Spot Colors

Cyan
0%

Magenta
0%

Yellow
0%

Black
100%

PMS
2945C

Special Instructions:

Quality Requirement of Artwork and Quality
Clarification Process of Artwork Printing:
Meet Eng-QS-06&02

Helen of Troy
Creative Services
Marlborough, MA 01752 USA
+1 508 490 7000